

# CONTRAINDICATIONS

Any person falling into one or more of the categories listed below should not use the *IonCleanse*:

- ❖ Wearer of a pacemaker or any other battery operated or electrical implant.
- ❖ Patients or clients who do not respond to a muscle test, or who test weak to both polarities, with their feet in the water and with the unit turned on.
- ❖ Those on heartbeat regulating medication.
- ❖ Organ transplant recipients.
- ❖ Organ removed
- ❖ Anyone on medication, the absence of which would mentally or physically incapacitate them, such as psychotic episodes, seizures, etc.
- ❖ Patients who are suffering from pathologies and are under the allopathic care system.
- ❖ Patients suffering from pathologies and are under the allopathic care system

Furthermore, the following recommendations should be strongly considered:

- ❖ Many medications require that a blood level be maintained in order to be effective, for example, blood pressure medication. Try to schedule an *IonCleanse* session just before the taking of a medication so that the patient or client can maintain proper blood levels of their medication.
- ❖ Patients with low blood sugar should eat before receiving a session. The *IonCleanse* tends to lower blood sugar in diabetics and may do so with those who are hypoglycemic.
- ❖ It is recommended that pregnant women and breast-feeding mothers not bathe because the *IonCleanse* will pull toxins out of tissues and fat that **could** impair fetal development and **will** become part of the milk that is fed to a newborn infant.  
**Note:** Pregnant and lactating women have used the *IonCleanse* without harm, but we recommend being on the safe side and not bathing them.
- ❖ As a general rule, *IonCleanse* can be used with patients on dialysis, taking insulin, or experiencing congestive heart failure. Gentle detoxification will help the body to eliminate the toxins that the kidney and heart cannot eliminate on their own and will not interfere with medications or deplete insulin levels.
- ❖ Patients who have a metal joint implant can be bathed, but some have found the exposure to an electromagnetic field to be uncomfortable. In the event of discomfort, notify the doctor and the session will be stopped immediately.

I have read and fully understand the above contraindications

\_\_\_\_\_  
Signature

\_\_\_\_\_  
Date

\_\_\_\_\_  
Print Name